

CLEANING GUIDE FOR

HCL® Easy Read Vials Stand, Grippers & Sets (# 18837, 18838, 18839, 18835B, 18836B, 18900, 18832, 18833 & 18834)

Materials:

Tray: Polystyrene

Grippers: GS-30HTA Clear Silicone Base

Recommended Disinfectant Options:

19118 Sterile Isopropyl Alcohol, 16 oz.

19185 Sterile HYPO-CHLOR® 0.52% Wipes, 12 x 12

Alternate Disinfectant Options:

70% Isopropyl Alcohol

0.52% Hydrogen Peroxide

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.