

CLEANING GUIDE FOR HCL® Hazardous Drug Tray, 11 x 12 (#20480)

Materials:

Polystyrene

Recommended Disinfectant Options:

19118 Sterile Isopropyl Alcohol, 16 oz.

19174 Sterile ALCOH-WIPE® 70% Individual Wipes, 12 x 12

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7, Canister

19185 Sterile HYPO-CHLOR® 0.52% Wipes, 12 x 12

Alternate Disinfectant Options:

70% Isopropyl Alcohol

0.52% Hydrogen Peroxide

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com

© Health Care Logistics®, Inc. 2026