

CLEANING GUIDE FOR

HCL® Unit Dose Bins (#5300, 5305 & 5310)

Materials:

Polycarbonate

Recommended Disinfectant Options:

20142 Sterile HCL® Presaturated Polypropylene Wipes, 70% IPA, 9 x 11

17880 SatPax® Presaturated Cleanroom Wipes, Non-Sterile, 9 x 9

Alternate Disinfectant Options:

70% Isopropyl Alcohol Based Cleaners

3% Hydrogen Peroxide Based Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.