

CLEANING GUIDE FOR

EZ Grid® Applicator Frame (#6520)

Materials:

6061 Aluminum

Recommended Disinfectant Options:

19118 Sterile Isopropyl Alcohol, 16 oz.

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12

Alternate Disinfectant Options:

70% Isopropyl Alcohol

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.