

INSTRUCTIONS FOR HCL® Tamper Tuf™ Vials w/ Caps (#7741_7743_7748_7749)



Step 1. Securely hold the Tamper Tuf Vial, then push the thumb tab outward.



Step 2. Carefully pull back on the tear-away strip to avoid popping off the top. Hint: To start the tear, pull toward the bottom of the container and twist.



Step 3. Continue pulling the tear-away strip until it reaches the end.



Step 4. Remove the cap by grasping the ribbed edges, twisting and lifting while securely holding the vial with your other hand.