

CLEANING GUIDE FOR

HCL® Personal Tablet Crusher w/ Cap (#7335, 7338, 21294 & 21295)

Materials:

Polycarbonate & Polypropylene

Recommended Disinfectant Options:

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12

20142 Sterile HCL® Presaturated Polypropylene Wipes, 70% IPA, 9 x 11

Alternate Disinfectant Options:

70% Isopropyl Alcohol based cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.