

CLEANING GUIDE FOR

HCL® Mobile Tables (#19102, 19103, 19102, 20880, 21091 & 21092)

Materials:

Stainless Steel and Polypropylene, Thermoplastic rubber(casters)

Recommended Disinfectant Options:

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12

19118 Sterile Isopropyl Alcohol, 16 oz.

Alternate Disinfectant Options:

70% Isopropyl Alcohol

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com