

# **CLEANING GUIDE FOR HCL® Vial Trays, Plastic (#20793 & 20794)**

**Materials:**

HDPE

**Recommended Disinfectant Options:**

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 9, Canister

**Alternate Disinfectant Options:**

70 % Isopropyl Alcohol

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

**Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923**

GoHCL.com • hcl@GoHCL.com

© Health Care Logistics®, Inc. 2026