

CLEANING GUIDE FOR

HCL® Medication Prep Tray, Stainless Steel (#20656)

Materials:

Stainless Steel

Recommended Disinfectant Options:

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPS, 9 x 12

19118 Sterile Isopropyl Alcohol, 16 oz.

Alternate Disinfectant Options:

70% Isopropyl Alcohol

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com

© Health Care Logistics®, Inc. 2026