

HCL® Bottle/ Vial #	8010-10 HCL® Adapter Caps for liquids, 18mm/A	7884-10 HCL® Adapter Caps for Liquids, 22mm/C	8018-10 HCL® Adapter Caps for Liquids, 28mm/K	8014-10 HCL® Adapter Caps for Liquids, 33mm/G	8012-10 Adapter cap 20mm/B	7886-10 Adapter Cap 24mm/D	8019-10 Adapter Cap 24mm/M	7883-10 Adapter Cap 28mm/E	8013-10 Adapter Cap 28mm/F
19489	X	X	X	X	X	✓	✓	X	X
19490	X	X	X	X	X	✓	✓	X	X
19491	X	X	X	X	X	✓	✓	X	X
21074	X	X	✓	X	X	X	X	✓	✓
20833	X	X	✓	X	X	X	X	✓	✓
21316	X	X	X	X	X	✓	✓	X	X
21317	X	X	X	X	X	✓	✓	X	X
10289	X	X	✓	X	X	X	X	✓	✓
10290	X	X	✓	X	X	X	X	✓	✓
10291	X	X	✓	X	X	X	X	✓	✓
10292	X	X	✓	X	X	X	X	✓	✓
10293	X	X	✓	X	X	X	X	✓	✓
10078	✓	X	X	X	X	X	X	X	X
10079	✓	X	X	X	X	X	X	X	X
10080	X	✓	X	X	X	X	X	X	X

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	8010-10 HCL® Adapter Caps for liquids, 18mm/A	7884-10 HCL® Adapter Caps for Liquids, 22mm/C	8018-10 HCL® Adapter Caps for Liquids, 28mm/K	8014-10 HCL® Adapter Caps for Liquids, 33mm/G	8012-10 Adapter cap 20mm/B	7886-10 Adapter Cap 24mm/D	8019-10 Adapter Cap 24mm/M	7883-10 Adapter Cap 28mm/E	8013-10 Adapter Cap 28mm/F
19368	✗	✗	✗	✗	✗	✓	✓	✗	✗
19369	✗	✗	✗	✗	✗	✗	✗	✗	✗
19365	✓	✗	✗	✗	✗	✗	✗	✗	✗
19366	✗	✗	✗	✗	✓	✗	✗	✗	✗
19367	✗	✗	✗	✗	✓	✗	✗	✗	✗
19370	✓	✗	✗	✗	✗	✗	✗	✗	✗
19450	✓	✗	✗	✗	✗	✗	✗	✗	✗
19451	✓	✗	✗	✗	✗	✗	✗	✗	✗
19452	✓	✗	✗	✗	✗	✗	✗	✗	✗
19453	✓	✗	✗	✗	✗	✗	✗	✗	✗
19454	✓	✗	✗	✗	✗	✗	✗	✗	✗
19455	✓	✗	✗	✗	✗	✗	✗	✗	✗
10407	✗	✗	✓	✗	✗	✗	✗	✓	✓
10408	✗	✗	✓	✗	✗	✗	✗	✓	✓
10410	✗	✗	✓	✗	✗	✗	✗	✓	✓
10409	✗	✗	✓	✗	✗	✗	✗	✓	✓

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	8010-10 HCL® Adapter Caps for liquids, 18mm/A	7884-10 HCL® Adapter Caps for Liquids, 22mm/C	8018-10 HCL® Adapter Caps for Liquids, 28mm/K	8014-10 HCL® Adapter Caps for Liquids, 33mm/G	8012-10 HCL® Adapter Cap 20mm/B	7886-10 HCL® Adapter Cap 24mm/D	8019-10 HCL® Adapter Cap 24mm/M	7883-10 HCL® Adapter Cap 28mm/E	8013-10 HCL® Adapter Cap 28mm/F
5703-04	X	X	X	X	✓	X	X	X	X
17868-01	X	X	X	X	✓	X	X	X	X
17869-01	X	X	X	X	✓	X	X	X	X
21314	X	X	X	X	✓	X	X	X	X
21315	X	X	X	X	✓	X	X	X	X
18580A	X	X	X	✓	X	X	X	X	X
18580G	X	X	X	✓	X	X	X	X	X
18582A	X	X	X	✓	X	X	X	X	X
18582G	X	X	X	✓	X	X	X	X	X

26-326 - ©Health Care Logistics®, Inc. 2026 REV-02



= Fits

HCL® Bottle/ Vial #	19393 Press in Bottle Adapter 24mm	19394 Press in Bottle Adapter 28mm	17282 Press In Bottle Adapter 20mm	17283 Press In Bottle Adapter 24mm	17284 Press In Bottle Adapter 28mm	7849 (31) HCL® Tapered Bottle Adapter	21142 NeoMed® Pharmacy Caps, 18mm	21143 NeoMed® Pharmacy Caps, 20mm	21144 NeoMed® Pharmacy Caps, 22mm	21145 NeoMed® Pharmacy Caps, 24mm	21146 NeoMed® Pharmacy Caps, 28mm
19489	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
19490	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
19491	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
21074	✗	✓	✗	✗	✓	✓	✗	✗	✗	✗	✓
20833	✗	✓	✗	✗	✓	✓	✗	✗	✗	✗	✓
21316	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
21317	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
10289	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10290	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10291	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10292	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10293	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10078	✓	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗
10079	✓	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗
10080	✗	✗	✓	✗	✗	✗	✗	✗	✓	✗	✗

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	19393 Press in Bottle Adapter 24mm	19394 Press in Bottle Adapter 28mm	17282 Press In Bottle Adapter 20mm	17283 Press In Bottle Adapter 24mm	17284 Press In Bottle Adapter 28mm	7849 (31) HCL® Tapered Bottle Adapter	21142 NeoMed® Pharmacy Caps, 18mm	21143 NeoMed® Pharmacy Caps, 20mm	21144 NeoMed® Pharmacy Caps, 22mm	21145 NeoMed® Pharmacy Caps, 24mm	21146 NeoMed® Pharmacy Caps, 28mm
19368	✓	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗
19369	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19365	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19366	✗	✗	✓	✗	✗	✗	✗	✓	✗	✗	✗
19367	✗	✗	✓	✗	✗	✗	✗	✓	✗	✗	✗
19370	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19450	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
19451	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
19452	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
19453	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
19454	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
19455	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
10407	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10408	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10410	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓
10409	✗	✗	✗	✗	✓	✓	✗	✗	✗	✗	✓

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	19393 Press in Bottle Adapter 24mm	19394 Press in Bottle Adapter 28mm	17282 Press In Bottle Adapter 20mm	17283 Press In Bottle Adapter 24mm	17284 Press In Bottle Adapter 28mm	7849 (31) HCL® Tapered Bottle Adapter	21142 NeoMed® Pharmacy Caps, 18mm	21143 NeoMed® Pharmacy Caps, 20mm	21144 NeoMed® Pharmacy Caps, 22mm	21145 NeoMed® Pharmacy Caps, 24mm	21146 NeoMed® Pharmacy Caps, 28mm
5703-04	X	X	✓	X	X	✓	X	X	✓	X	X
17868-01	X	X	✓	X	X	X	X	X	✓	X	X
17869-01	X	X	✓	X	X	X	X	X	✓	X	X
21314	X	X	✓	X	X	✓	X	✓	X	X	X
21315	X	X	✓	X	X	✓	X	✓	X	X	X
18580A	X	X	X	X	X	X	X	X	X	X	X
18580G	X	X	X	X	X	X	X	X	X	X	X
18582A	X	X	X	X	X	X	X	X	X	X	X
18582G	X	X	X	X	X	X	X	X	X	X	X

26-326 - Health Care Logistics®, Inc. 2026 REV-02



= Fits

HCL® Bottle/ Vial #	20227 (31) SealSafe™ Pierced Bottle Adapters 20mm	20228 (31) SealSafe™ Slit Bottle Adapters 20mm	20229 (31) SealSafe™ Slit Bottle Adapters 24mm	20230 (31) SealSafe™ Slit Bottle Adapters 28mm	17727 (31) SealSafe™ Pierced Bottle Adapters 24mm	18761 (31) SealSafe™ Pierced Bottle Adapters 28mm	20312 ENFit® Bottle Adapter 12-15mm	20313 ENFit® Bottle Adapter 14-16.5mm	20314 ENFit® Bottle Adapter 17-20.5mm	20315 ENFit® Bottle Adapter 19-23mm	21052 ISOSAF™ ENFit®/Oral Stepped Bottle Adapt.
19489	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗	✓
19490	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗	✓
19491	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗	✓
21074	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
20833	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
21316	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗	✓
21317	✗	✗	✓	✗	✓	✗	✗	✗	✓	✗	✓
10289	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10290	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10291	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10292	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10293	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10078	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
10079	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗
10080	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	20227 (31) SealSafe Pierced Bottle Adapters 20mm	20228 (31) SealSafe Slit Bottle Adapters 20mm	20229 (31) SealSafe™ Slit Bottle Adapters 24mm	20230 (31) SealSafe™ Slit Bottle Adapters 28mm	17727 (31) SealSafe™ Pierced Bottle Adapters 24mm	18761 (31) SealSafe™ Pierced Bottle Adapters 28mm	20312 ENFit® Bottle Adapter 12-15mm	20313 ENFit® Bottle Adapter 14-16.5mm	20314 ENFit® Bottle Adapter 17-20.5mm	20315 ENFit® Bottle Adapter 19-23mm	21052 ISOSAF™ ENFit/Oral Stepped Bottle Adapt.
19368	✗	✗	✓	✗	✓	✗	✗	✗	✗	✗	✓
19369	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19365	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19366	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✗
19367	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✗
19370	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19450	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19451	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19452	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19453	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19454	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
19455	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
10407	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10408	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10410	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓
10409	✗	✗	✗	✓	✗	✓	✗	✗	✗	✓	✓

26-326 - Health Care Logistics®, Inc. 2026 REV-402



= Fits

HCL® Bottle/ Vial #	20227 (31) SealSafe™ Pierced Bottle Adapters 20mm	20228 (31) SealSafe™ Slit Bottle Adapters 20mm	20229 (31) SealSafe™ Slit Bottle Adapters 24mm	20230 (31) SealSafe™ Slit Bottle Adapters 28mm	17727 (31) SealSafe™ Pierced Bottle Adapters 24mm	18761 (31) SealSafe™ Pierced Bottle Adapters 28mm	20312 ENFit® Bottle Adapter 12-15mm	20313 ENFit® Bottle Adapter 14-16.5mm	20314 ENFit® Bottle Adapter 17-20.5mm	20315 ENFit® Bottle Adapter 19-23mm	21052 ISOSAF™ ENFit/Oral Stepped Bottle Adapt.
5703-04	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✗
17868-01	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✗
17869-01	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✗
21314	✓	✓	✗	✗	✗	✗	✓	✗	✗	✗	✓
21315	✓	✓	✗	✗	✗	✗	✗	✓	✗	✗	✗
18580A	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
18580G	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
18582A	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓
18582G	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓

The following adapters don't fit any of the HCL items
 19392 Press in Bottle Adapter 22mm
 17285 Press In Bottle Adapter 33mm

19541 Adapter Cap 38mm/
 20316 ENFit® Bottle Adapter 24-26mm
 21145 NeoMed® Pharmacy Caps, 24mm

21147 NeoMed® Pharmacy Caps, 33mm
 21148 NeoMed® Pharmacy Caps, 38mm

