

CLEANING GUIDE FOR

Warming Cabinets (#4833, 4834, 4857 & 4858)

Materials:

Laminated Partical Board

Recommended Disinfectant Options:

pH-neutral Cleaner

Alternate Disinfectant Options:

Diluted Vinegar and Water Solution

1. Use a 1:4 mixture (1 part vinegar to 4 parts water) and apply it with a well-wrung microfiber cloth to avoid excess moisture, which can damage the surface.

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.